

## SURGICAL KITS

### IMPLANCE SURGICAL KITS

BONE - TISSUE LEVEL



### SHORT IMPLANT SURGICAL KITS



### AGGRESSOR SURGICAL KITS



FDA CE  
1984



AGS Medikal Ürünleri İth. İhr. Tic. A.Ş.  
Katip Salih Sok. No: 111 34718 Koşuyolu / İstanbul / Türkiye  
T +90 216 546 00 78 F +90 216 546 00 79

implance.com

July 2022



IMPLANCE  
Dental Implant System



## Implance Dental Implant System

**Implance** Dental Implant System, has played a very important role in the growth of the dental implantology sector in Turkey by developing new products of superior quality tailored to the needs of the industry.

**Implance**, with its investments in technological production and qualified staff, can be supplied in Turkey and all over the world in a short time, it provides the opportunity to find the best quality-price ratio.


**Implance**, has achieved world-wide recognition with many successful works from the moment of its establishment and became the first **FDA\*** approved dental implant brand in Turkey.

Approved by the **FDA**, **Implance** has become one of the most trusted brands in the world by entering the US market with its outstanding product quality and special service network.


\* American Food and Drug Administration




## BONE LEVEL IMPLANT




Length (mm)	Code
8.0	BMFSR3308
10.0	BMFSR3310
12.0	BMFSR3312
14.0	BMFSR3314




Length (mm)	Code
8.0	BFHR3708
10.0	BFHR3710
12.0	BFHR3712
14.0	BFHR3714



Length (mm)	Code
8.0	BFHR4308
10.0	BFHR4310
12.0	BFHR4312
14.0	BFHR4314




Length (mm)	Code
8.0	BFHR4808
10.0	BFHR4810
12.0	BFHR4812
14.0	BFHR4814




Length (mm)	Code
8.0	BWFHR5508
10.0	BWFHR5510
12.0	BWFHR5512


## TISSUE LEVEL IMPLANT




Length (mm)	Code
8.0	TFHR3708
10.0	TFHR3710
12.0	TFHR3712
14.0	TFHR3714



Length (mm)	Code
8.0	TFHR4308
10.0	TFHR4310
12.0	TFHR4312
14.0	TFHR4314



Length (mm)	Code
8.0	TFHR4808
10.0	TFHR4810
12.0	TFHR4812
14.0	TFHR4814




Length (mm)	Code
8.0	TWFHR5508
10.0	TWFHR5510
12.0	TWFHR5512


## SHORT IMPLANT

In cases where bone level is insufficient in the upper and lower jaws, and when neither the mandibular nerve nor the sinus allow for regular implants, short implants may be an alternative solution instead of resorting to advanced surgery techniques.


### BONE LEVEL



Length (mm)	Code
5.0	BMFSR4305
6.0	BMFSR4306




Length (mm)	Code
5.0	BMFSR4805
6.0	BMFSR4806




Length (mm)	Code
5.0	BMFSR5505
6.0	BMFSR5506


### TISSUE LEVEL



Length (mm)	Code
4.0	TFSR4304
5.0	TFSR4305
6.0	TFSR4306



Length (mm)	Code
4.0	TFSR4804
5.0	TFSR4805
6.0	TFSR4806




Length (mm)	Code
4.0	TWFSR5504
5.0	TWFSR5505
6.0	TWFSR5506


## AGGRESSOR IMMEDIATE LOADING

Aggressor implant offers the advantage of immediate loading. It does not cause compression at the neck area due to its aggressive thread design and sharper threads towards the apex.


Aggressor implant body design with increased taper towards its apex, provides high primary stability in the extraction sockets and D3-D4 bone types.



Length (mm)	Code
10.0	BFAR3710
12.0	BFAR3712
14.0	BFAR3714



Length (mm)	Code
10.0	BFAR4310
12.0	BFAR4312
14.0	BFAR4314



Length (mm)	Code
10.0	BFAR4810
12.0	BFAR4812
14.0	BFAR4814